

CideTest 1 gt™

Glutaraldehyde Potency Test Strips | 1.0–2.5%

Assurance your products are safe and effective

- Glutaraldehyde test strips are necessary to assess the concentration of glutaraldehyde in high level disinfecting and sterilizing solutions. Using test strips verify the percentage of glutaraldehyde in the solution.
- Quality control tool monitors the effectiveness of Germiphene's Gluterate (2.0% glutaraldehyde) solution during reprocessing of medical equipment and instruments.
- High-Level Disinfectants and Sterilants must maintain the minimum effective concentration (MEC) level of 1.5% glutaraldehyde or greater during the 28-day use life to effectively sterilize and high-level disinfect.



Use with Gluterate!

2% Glutaraldehyde Solution




You should test the solution on the basis of its use. For instance, if you use the solution daily, it must be tested daily and the results recorded in your high level disinfectant/sterilization log.

Easy & Quick to Use

Dip strip in test solution for 1 second, or apply a drop of test solution onto the reagent pad. Remove excess test solution by gently shaking the strip or by touching the side of the strip in a paper sorbent. At 30 seconds, compare the strip colour with the colour blocks.

Fast & Clear Results

Glutaraldehyde Concentration



For more information please contact your District Sales Manager or call Germiphene at 1.800.265.9931 | www.germiphene.com

R0