

Disinfectant and cleaner for use on non-invasive devices, equipment and surfaces in medical, surgical and dental health care facilities such as stethoscopes, stainless steel or plastic trays and alligator clips.

GamaCide effectively kills the following microorganisms at room temperature (20 C) with a 1 minute contact time: *Mycobacterium terrae* ATCC 15755, *Staphylococcus aureus* ATCC 6538, *Salmonella choleraesuis* ATCC 10798, *Pseudomonas aeruginosa* ATCC 15442, *Escherichia coli* NCTC 8196, *Enterococcus hirae* ATCC 10541, Poliovirus 1 strain Chat ATCC VR-1562 and Influenza A viruses. Effectively inactivates SARS-Cov-2 (COVID 19), Herpes Simplex 1 virus, Human immunodeficiency virus type 1, Human hepatitis B virus in 2 minutes. A 3-minute contact time is required to kill *Trichophyton mentagrophytes* ATCC9533.

**DIRECTIONS:** Appropriate glove protection is recommended prior to use of product. Use where control of the hazards of cross-contamination is of concern. Remove large cover. Remove safety seal. Locate first wipe in center of the roll and thread through eyelet of cover. Replace cover. Pull each wipe out at an angle from the cap. Close cap after each use. **DO NOT SWALLOW.**

**GENERAL DISINFECTION:** Thoroughly clean the surface of heavy and organic soil prior to disinfection. Wipe the pre-cleaned area to be disinfected. Allow surface to remain wet for 3 minutes for broad spectrum disinfection. Surgical instruments must be sterilized according to manufacturer's instructions after cleaning with GamaCide<sup>3</sup> wipes.

**PERSONAL PROTECTION:** Wear appropriate barrier protection such as gloves, gowns, masks or eye coverings. Cleaning Procedure: Blood and other body fluids must be thoroughly cleaned from surfaces/objects before disinfection. Wipe the pre-cleaned area to be disinfected. Allow surface to remain wet for 3 minutes for broad spectrum disinfection.

**INFECTIOUS MATERIALS DISPOSAL:** Cleaning materials used that may contain blood/body fluids should be autoclaved and/or disposed of in accordance with local regulations for infectious materials disposal.

**CONTRAINDICATIONS:** Not for use in hospital nurseries or on children's toys. Professional use only. Do not use to disinfect invasive devices. Invasive devices must be disinfected or sterilized prior to use.



**Disinfectant for Non-Invasive Devices, Equipment, and Surfaces**  
**Bactericidal • Fungicidal • Virucidal • Tuberculocidal**  
**Medical/Surgical/Dental Disinfectant Cleaner**  
For Professional Use Only. Ready to Use.

**ACTIVE INGREDIENTS:**  
Benzethonium Chloride.....0.35% w/w  
Isopropyl Alcohol.....20.00% w/w

**160**  
15x17cm  
Premoistened  
Cloths

MD REF EU-GAMA-SW

UDI



(01) 0 0850046 43639 0  
(10) LOT123  
(17) 260309

LOT ???????

USE BY: YYYY-MM-DD

**STORAGE AND DISPOSAL:** Do not contaminate food or feed by storage or disposal. Storage: Store in cool, well ventilated area. Product Disposal: Dispose of wipe in trash after use. Unused contents to be disposed of in accordance with local/regional/national regulations. Container Disposal: Triple rinse empty container, then offer for recycling or reconditioning; or puncture and dispose of in a sanitary landfill.

**FLAMMABLE:** Do not use in presence of open flame or spark.  
**CAUTION: KEEP OUT OF REACH OF CHILDREN.**

Warning! Causes eye and skin irritation. Avoid contact with eyes, skin, or clothing. Wear protective gloves, clothing, goggles, face protection. Wash thoroughly with soap and water after handling. Remove contaminated clothing and wash clothing before reuse.

**FIRST AID:** For more information, refer to SDS. **IF ON SKIN:** Take off contaminated clothing. Rinse skin immediately with water. **IF IN EYES:** Rinse eyes with water for several minutes. Remove contacts, if worn, and continue rinsing until advised to stop by the Poisons Information Centre or a doctor, or for at least 15 minutes. If skin irritation occurs or eye irritation persists, get medical attention. Have the product container/label with you when calling a Poison Control Centre or doctor, or going for treatment.

Young Microbrush Ireland Ltd.  
Clogherane, Dungarvan  
Co. Waterford, X35 VE02, Ireland

EU-GM3-LB-SW R1



0050

EC REP CEpartner4U  
Esdoornlaan 13, 3951 DB Maam, Netherlands

Manufactured by: **Germiphene Corporation**  
1379 Colborne St. E. | Brantford ON | Canada N3T 5M1  
Toll Free Phone 1.800.265.9931

For a list of additional translations, please go to:  
<https://germiphene.com/product-resources/>

MADE IN CANADA

13 in

3.24 in