



Dental Waterline Treatment Maintenance & Shock

Specifically designed for use in independent water bottle systems.

A pump a day keeps the bacteria and biofilm away

- Biofilm & bacteria control
- No harsh chemicals
- Odourless, colourless, tasteless, pH neutral
- Non foaming, consistent results every time
- Tested on various waterline tubing with no negative effects or corrosion
- One E2O™ bottle lasts approximately 4–6 months per dental chair

Eco-friendly solution, E2O™ breaks down into water and oxygen
Non-toxic to equipment or patients

VS

- No waiting for tablets to dissolve, no residue/clogging from undissolved tablets
- No patient complaints about bad taste or grit from tablet residue

Take charge

of your office water today

Easy dose pump.
No measuring, no mess.

Convenience at an economical price

\$SAVE:

Cost per use

\$0.51*

REGULAR TREATMENT

700 mL
2L


SHOCK TREATMENT



Water sample taken from an ultrasonic scaler



Before E2O™



After E2O™

▲ Tests performed showed a significant reduction of bacteria in water sources after the use of E2O™.

SAFE to use

A whitening strip, directly placed in your mouth, contains **500X** more hydrogen peroxide than E2O™ in your waterlines!*

Interested in waterline testing?

We ♥ to educate! Learn about dental waterlines – the challenges and solutions. Visit germiphene.com for our complimentary on demand webinars or book an office session with your rep today.

E2O 250mL.....E2O-M

Powered by
Huwa-San
PEROXIDE TECHNOLOGY

*Based on using the daily treatment dosage of E2O™ in your waterline unit.
* Cost per use based on 700 mL water bottle and daily treatment dosage.



1.800.265.9931

www.germiphene.com



Which testing is best for your office?

introducing our new online waterline app

GET MORE WITH E|One™ + E|Line™

Schedule waterline tests + record results with ease and accuracy.

Ask us how to get started today!



Water Quality Check Kit

2 Test vials:


- Total Bacteria + Fungi
Test detection limit 200 cfu/mL
- *E. coli* + *faecal coliform*
Test detection limit 1 cfu/mL

IN HOUSE TESTING




RESULTS
in 24 hours
using a 37°C Dry
Block Incubator

Easy to read,
colour coded results



E|Line™ Complete Control


- Complete control, all in one box



18 tests

SKU: E2O-KIT


E|Line™ Water Quality Check Kit



18 tests

SKU: ELINE-W-KIT

E|Line™ Water Quality Check



12 tests

SKU: ELINE-VX2



Water Quality Indicator

- Broad detection
- One vial tests for bacteria, biofilm & fungi

Test detection limit 500 cfu/mL



FAST RESULTS

Convenient in house testing



RESULTS
quick read at
15 minutes



RESULTS
final read at
1 hour

EASY TO READ

Clear results



NO

Incubator required



Clean water



Contaminated
water



Highly
contaminated water

12 TEST KITS PER BOX

SKU: E1-12X20ML-SYR-KIT


Each Test Kit Includes:



Disposable Filter



Saline Water Vial



Disposable Syringe



E|One Test Vial