

SURGICAL MILK

Concentrated Instrument Lubricator and Rust Inhibitor

Get the Most From Your Instruments!

- Keeps instruments functional and prevents dulling
- Clinically tested to be steam penetrable and steam sterilizable
- Effective in steam, dry or chemical vapour sterilizers
- Reduces instrument repair and replacement costs



It's Effective	Bacteriostatic.
	Prevents spotting, staining and instrument rust. Keeps instruments shiny and prevents dulling.
	Effective in steam, dry or chemical vapour sterilizers.
	Routine use reduces cost of instrument repair and replacement.
	Clinically tested to be steam penetrable and steam sterilizable. Will not interfere with the microbial killing efficacy of autoclaving.
It's Convenient	Not sticky, oily or greasy.
	1 part solution to 5 parts water can be used for up to 7 days.
	Can be used on a variety of metals: steel, carbide, tungsten, stainless steel and aluminum.



2L

Item #:
SM-S 2L



4L

Item #:
SM-S 4L



made in
Canada

4L Bottle of Concentrate = 20L of Ready-to-Use Solution



Use During Autoclaving on

- Dental/Surgical Instruments
- Hinges
- Box Locks
- General Equipment

To order Surgical Milk or for more information, contact your Account Manager or Germiphene at 1.800.265.9931 | info@germiphene.com | www.germiphene.com

REV 01 072024

GERMIPHENE

A Division of Young Innovations