

GLEAM II

Ultrasonic Cleaning Solution

A pleasantly scented cleaning solution that will **make your instruments shine.**

- Bacteriostatic
- Pleasantly scented
- Non-ammoniated formulation
- Will not harm stainless steel instruments
- Concentrated solution - 100 mL makes 1 L of use solution



Convenient	<p>Concentrated with a dilution rate of 1 Part Gleam II to 10 parts water saves on storage space</p> <p>Easy to use 2L tip and measure bottle is available which prevents any waste.</p>
Effective	<p>Specifically developed for optimal cleaning efficiency in all types of ultrasonic cleaning units</p> <p>Powerful detergents and surfactants that will clean dental instruments quickly and safely by removing blood and organic soils of all types</p>
Safe	<p>Will not harm stainless steel instruments</p> <p>Unique dispensing system avoids any spills</p>



2 L

Item #: GLE



10 L

Item #: GLE-TP



made in
Canada

INSTRUMENT REPROCESSING BACKED UP?

ET. FOAM | Enzymatic Soaking Foam

- Keeps soiled instruments moist
- Ideal for use prior to instrument/washer ultrasonic cleaning units
- Enzyme boost

To order Gleam II or for more information, contact your Account Manager or Germiphene at 1.800.265.9931 info@germiphene.com | www.germiphene.com



REV 01 062024

GERMIPHENE

A Division of Young Innovations