

# GLEAM II



## Ultrasonic Cleaning Solution

A pleasantly scented cleaning solution that will make your instruments shine.



Gleam II is a concentrated, non-ammoniated ultrasonic cleaning solution specially developed for optimal cleaning efficiency. Gleam II is designed to work in all types of ultrasonic cleaning units and is packaged in a convenient Tip 'n' Measure bottle for easy dispensing.

Gleam II's powerful detergents thoroughly clean dental instruments, quickly and safely removing blood, tissue, mucous and other organic soils. A unique blend of surfactants efficiently cleans away lipids, starches and inorganic debris and prevents them from re-depositing. The powerful formulation of Gleam II contains chelating agents to help restore instrument's original shine.

Gleam II contains special corrosion inhibitors and when used as directed will not harm stainless steel, hard plastics, rubber or glass. Gleam II leaves no residue on instruments or the ultrasonic cleaning tank.

- Bacteriostatic
- Pleasantly scented
- Non-ammoniated formulation
- Will not harm stainless steel instruments
- Concentrated solution - 100 mL makes 1 L of use solution

*Make your instruments shine!*



To order this or any Germiphene product, contact us by calling 1-800-265-9931 or via e-mail: [info@germiphene.com](mailto:info@germiphene.com)





# PRODUCT PROFILE

## GLEAM II

<b>Product</b>	Ultrasonic Cleaning Solution
<b>Description</b>	A concentrated, non-ammoniated formulation specifically developed for optimal cleaning efficiency in all types of ultrasonic cleaning units.
<b>Function</b>	<b>GLEAM II</b> is formulated to remove proteinaceous and other organic matter, such as blood, mucus, calculous, etc. from dental instruments.
<b>Packaging</b>	2 litre Tip 'n' Measure bottle and 10 litre container.

Features	Benefits
<i>Easy to use</i>	Simply tip the bottle to fill the upper chamber with the desired amount and pour into ultrasonic tank. This unique dispensing system prevents wastes and spills.
<i>Concentrate</i>	Dilute 1 part of concentrate with 10 parts water.
<i>Superior Cleaning Action</i>	<b>GLEAM II</b> 's formula with powerful detergents and surfactants will clean dental instruments quickly and safely by removing blood and organic soils of all types.
<i>Non-Corrosive</i>	Will not harm stainless steel instruments
<i>Fresh Scent</i>	<b>GLEAM II</b> has a pleasant scent
<i>Reduces Risk Of Cross-Contamination</i>	Ultrasonic cleaning is a critical step in the total infection control process. Prior to sterilization, critical instruments should always be cleaned ultrasonically (to remove all debris and contaminants), thus reducing the risk of cross-contamination associated with manual cleaning.
<i>Protection and Peace of Mind</i>	<b>GLEAM II</b> Ultrasonic Solution plays a vital role in protecting dental professionals by providing superior cleaning results.

If you have any questions about **GLEAM II**, contact your local area sales representative or call **GERMIPHENE** Corporation toll-free at 1.800.265.9931  
Website: [www.germiphene.com](http://www.germiphene.com) E-mail: [admin@germiphene.com](mailto:admin@germiphene.com)

