


## 2% Gluteraldehyde solution

- · Environmentally friendly
- One-step solution for Glutaraldehyde and OPA disposal
- Changes harmful environmental chemicals into non toxic, non hazardous waste
- · Shows deactivation by colour change



Effective	A proprietary organic compound that <b>neutralizes solution before disposal.</b>
	One pouch will neutralize 1L of 2%–2.5% Glutaraldehyde or 4 L of 0.55% OPA solution.
	Gluterate is reusable for up to 30 days.
	Gluterate contains a <b>highly effective surfactant system</b> that ensures complete wetting. It is also <b>extremely active in the presence of soaps and organic soil.</b>
Easy	Ready-to-use. Simply pour GLUTE NEUT packet into solution. No mixing, no waste.
	GLUTE NEUT changes Glutaraldehyde or OPA solutions from yellow to red when neutralized.
Safe	The natural oxidization process will neutralize solutions for easy disposal. GLUTE NEUT does not affect organisms used in sewage treatment plants and will not harm water systems or aquatic life.
	<b>GLUTE NEUT eliminates vapours caused by agitation of solution during disposal.</b> Use of GLUTE NEUT eliminates eye and skin irritations resulting from spills or splashes when disposing of solution.

GLUTE NEUT	COMPETITIVE LIQUIDS
Stackable package, minimal storage space	960mL (32 ounces) bottle.
1 Box (16 x 50 g packets) neutralizes 16L glutaraldehyde or 64 L OPA.	1 Bottle (960 mL) neutralizes 7.5 L glutaraldehyde or 60.5L OPA.
Powder	Liquid
Minimal risk of spills, splatter and mess	Higher risk of spills, splatter and mess vs powder
No measuring, no waste	Liquids must be measured accurately.





## **Try This Valuable Support Product**

## **GLUTERATE**

## **Sterilizing and Disinfecting Solution**

10 Minute TB disinfection, 10 minute high level disinfection.

Size:4L Bottle | Item# GLUT

To order GLUTE NEUT or for more information, contact your Account Manager or Germiphene at 1.800.265.9931 info@germiphene.com | www.germiphene.com

Rev 01 072024