

Gluterate

"GLUT"

Sterilizing and Disinfecting Solution Solution de stérilisation et de désinfection

- 10 Minute TB disinfection, 10 minute high level disinfection
- 10 Hour sterilization
- Bactericidal, Fungicidal, Tuberculocidal, Virucidal,
- Sporicidal
- Ready to use solution — No activation required
- Chemically stable —No degradation at high temperatures
- Fresh mint scent masks aldehyde odour



made in
Canada

Effective	A stable 2% glutaraldehyde solution that meets all requirements as a superior sterilizing/disinfecting solution.
	Provides complete sterilization or high-level disinfection of semi-critical dental/medical instruments and devices by immersion in a manual (covered tray) system.
	Gluterate is reusable for up to 30 days based on cide strip concentration results.
	Gluterate contains a highly effective surfactant system that ensures complete wetting. It is also extremely active in the presence of soaps and organic soil.
Broad Spectrum	Biocidal to vegetative organisms including Pseudomonas aeruginosa, Staphylococcus aureus, Salmonella choleraesuis, pathogenic fungi as represented by Trichophyton mentagrophytes, and lethal to Mycobacterium species such as Mycobacterium terrae.
	Demonstrated broad spectrum virucidal activity against the golden standard for virucidal testing, Poliovirus 1. Proven to inactivate viruses in 10 minutes.
Safe	Unique anti-corrosion formulation will not damage instruments or devices (e.g. lenses, cements, plastics, rubber, steel, aluminum, brass) when used as directed.
	Safe to use when proper Personal Protective Equipment is used.

Try These Related Products

CIDE STRIPS

Glutaraldehyde Potency Test Strips

Easily test the concentration of glutaraldehyde in high level disinfecting.

100 strips

Item#: GT-232-100



GLUTE NEUT

Neutralizer for Glutaraldehyde

Environmentally friendly, one-step solution for the safe disposal of glutaraldehyde.

Size: 16 x 50 g pouches/box

Item# GLUTE-NEUT



To order Gluterate or for more information, contact your Account Manager or Germiphene at 1.800.265.9931 | info@germiphene.com | www.germiphene.com

Rev 01 072024

GERMIPHENE

A Division of Young Innovations