

HemaSeal-G

Economical, multi purpose desensitizing solution

Unbeatable Value! *Half the price*
2X *the competition!*

Avoid patient discomfort and compromised restorations

Controls dentinal fluid by combining with proteins in the dentin tubule and producing a physiologic barrier. Seals tubuli in seconds, eliminating fluid shifting, the major cause of post operative and cervical sensitivity.



NPN 02246832 | 10mL

- ✓ Apply to all prepared teeth where dentin is involved
- ✓ Use to desensitize in all restorative situations
- ✓ Long lasting relief in seconds
- ✓ Zero effect on bond strength
- ✓ Protect the tooth/restoration surface from bacterial growth
- ✓ 0.5% fluoride to reduce dentin sensitivity
- ✓ Contains 35% Hema, 5% glutaraldehyde
- ✓ Hema enhances the ability of hydrophilic bonding agents to penetrate into the tubules

Convenient dropper dispenser



Versatile

Use as a desensitizer alone in cervical abrasions and root sensitivity or in conjunction with all restorative material including amalgams, crowns, inlays, onlays, veneers or composite resin restorations. Use to prevent post-operative and micro leakage sensitivity. Allows for complete seating of all indirect restorations.

Easy to Use

One step chair side procedure. Simply apply one coat and remove excess with a short blast of air. No mixing or curing, no repeat applications required.



For more information please contact your District Sales Manager or call Germiphene at 1.800.265.9931

R0