



Single Dose

## Ultrasonic Cleaning Solution

- Quickly removes organic matter (mucous, blood, saliva)
- Sparkling clean results
- Bacteriostatic — inhibits microbial growth
- Contains rust inhibitors
- Free rinsing
- Excellent for cleaning burrs
- Safe for all types of ultrasonic cleaning units

### Super Concentrated

- Concentrated solutions = less storage space and more responsible choices
- Dilute 30mL per 4L water



1 box = 192L  
SuperSonic solution

48 x 30mL

## Liquid vs Tablets

- SuperSonic is available in single dose sachets - no mess, no mixing
- Start cleaning immediately, don't wait for slow dissolving tablets

## High cleaning efficacy, easy to use & cost-effective

SuperSonic outperforms the competition!

Product	Solution Preparation (mins)†	Cleaning Time (mins)	Total Cleaning Time (mins)
SuperSonic	0	7.03	7.03
Maxizyme Tabs™	6.17	14.24	20.41
Clean & Simple™	6.14	9.27	15.41

Maxizyme Tabs™ is a registered trademark of Henry Schein, Inc. Clean & Simple™ is a registered trademark of Tuttnauer.

## Instrument reprocessing backed up?



### ET.FOAM | Enzymatic Transportation Foam

- Keeps soiled instruments moist
- Ideal for use prior to instrument washers/ultrasonic cleaning units
- Enzyme boost

REMEMBER to rinse instruments thoroughly before placing in the instrument washer/ultrasonic cleaning unit when using this type of product!



Why pay more? Competitors need at least twice as much time to work and cost more!

TOSI test performed. TOSI = soil load that represents coagulated blood dried on a metal test strip placed between glass plates. TOSI simulates instruments with difficult to clean parts like crevices and hinges.

†Time for tablets to fully dissolve.

To order SuperSonic or for more information, contact your District Sales Manager or Germiphene at 1.800.265.9931 | info@germiphene.com | www.germiphene.com

