




DON'T RISK CROSS- CONTAMINATION



- + Reduces cross-contamination of viruses and bacteria
- + Safer alternative to autoclavable tip
- + Instant Dry Air
- + Powerful, consistent airflow
- + Save time and money


Contact your Germiphene account manager or our Customer Service team at **1-800-265-9931** or **orders@germiphene.com**.

Visit our website at **www.germiphene.com**.


QUALITY YOU CAN DEPEND ON

Experience the proven design, superior performance and better hygiene of Crystal Tip®

Make the Switch Today!



Available in Assorted Rainbow



Beveled Edge
Patient comfort
and less wear and
tear on O-rings

Two-Piece Construction
Complete air/water separation

Quad Air Core
Powerful airflow and
superior atomization

Tip-Lock Groove
Secure Fit, preventing
accidental ejection



DISPOSABLE AIR/WATER SYRINGE TIPS

Assorted Rainbow

250ct	CT1200
1500ct	CT2100

