## GERMISAFE<sup>TM</sup>

## **Process Challenge Device**

## **CONSISTENT RESULTS YOU CAN TRUST!**

The Germisafe™ Process Challenge Device (PCD) saves time and eliminates variability associated with "test packs" that are manually prepared in-house.

Designed to simulate the most difficult instruments to sterilize, the Germisafe™ PCD provides a defined resistance to the sterilization process and works to ensure effective sterilization of ALL items in the load. The Germisafe™ PCD is validated to perform equivalent to ANSI/AAMI 16-towel test, the universally accepted gold standard for process challenge devices.\*



- ✓ The only validated, reusable PCD on the market today
- ✓ Simulates the most difficult instrument(s) to sterilize
- ✓ Intuitive flip-top design for easy loading and viewing of results
- ✓ Provides a consistent, reproducible challenge in every load
- ✓ For Gravity and Dynamic Air Removal Cycles

Germisafe™ Process Challenge Device......PCD-GERM

\*Tested per standards outlined in the American National Standards Institute (ANSI) and Association for the Advancement of Medical Instrumentation (AAMI), *ANSI/AAMI ST79: 2012 Comprehensive guide to sterilization and sterility assurance in health care facilities.* Compatible with the following indicators: Germisafe™ Biological Indicator, Bionova Ultra Rapid Biological Indicator, Steamline™ Type 5 Chemical Integrators, and Vaporline® Type 5 Chemical Integrators

Performance validated for 600 cycles (~6 months)

Rev. 1 / 03.2025



NEW

BE SURE TO ASK ABOUT OUR IPAC EDUCATION OPTIONS!

www.germiphene.com | 1.800.265.9931