

NPN 01910531

Topical 1.23% Fluoride Gel

75 Uses Available in Every Bottle.

- Rapid fluoride update in 60 seconds •
- Recommended Adult Fluoride treatment¹ .
- Use with a tray or paint on method •

Gluten Alcohol Dye Free Free Free

PATIENT CHALLENGES

• Caries Age 6+* Adult Coronal Caries

To arrest or reverse noncavitated carious lesions on occlusal, facial or lingual surfaces of permanent and primary teeth the ADA Expert panel suggests clinicians use 1.23% acidulated phosphate fluoride gel (application every 3-6 months).¹

	It Works	A topical Acidulated Phosphate Fluoron (APF) anti-caries preparation gel containing 1.23% fluoride ion from sodium fluoride and hydrofluoric acid with low pH		
		Recommended for children and adults without composite and porcelain dental restorations, crowns or bridges		STRAWBERRY Item #: G6-S made in Canada
		Adheres closely to the enamel and interproximal areas		
	lt's Patient Friendly	Thixotropic gel minimizes risk of patient gagging/ingestion, maximizes contact with all tooth surfaces		
		Patient pleasing flavours		

To order 60 Second Gel or for more information, contact your Account Manager or Germiphene at 1.800.265.9931 | info@germiphene.com | www.germiphene.com

1.ADA Evidence-Based Clinical Practice Guideline on Nonrestorative Treatments for Carious Lesions: A report from the American Dental Association, 2018 * Exercise caution with younger children. Fluoride gels should not be swallowed.

REV 01 072024

450mL

MINT

Item #: G6-SP

GERMIPHENE

A Division of Young Innovations