

# TRÉS CLEAN

## Concentrated Impression Tray & Instrument Cleaner

- Loosens any remaining traces of alginates and zinc-oxide eugenol
- Mild cleaning formula leaves trays and instruments sparkling clean
- Non rusting formula is safe for all non corroding metal impression trays and instruments
- No messy powder to measure and dissolve, no dust



FRESH  
LEMON  
SCENT

1L

Cleans significantly better than the competition!


% Cleaned After 30 Minute Soak\*



Cleans better and **5x FASTER**  
than the old formula



COMPETITOR A  
(POWDER)



COMPETITOR B  
(POWDER)



COMPETITOR D  
(POWDER)

Clean  
responsibly

PHOSPHATE  
FREE

clean responsibly

Liquid  
formula!

1L of  
Trés Clean  
makes  
**20L**  
of solution



**FAST  
ACTING**



\*Data on file

1.800.265.9931 | [www.germiphene.com](http://www.germiphene.com)

